


Operating Table EXAFLEX 3000

Technical Data

Dimensions



Max. patient weight	160 kg
Weight (incl. head and leg plate)	165 kg
Height max.	1070 mm
Height min.	770 mm
Total width (incl. accessories rails)	590 mm
Width of bed surface	540 mm
Length	total 2000 mm
	Back plate 620 mm
	Seat plate 425 mm
	Leg plate 620 mm

Adjustment

Stroke	300 mm
Back plate, raise	75°
Back plate, lower	45°
Leg plate, raise	20°
Leg plate, lower	90°
Trendelenburg	30°
Reverse-Trendelenburg	30°
Lateral tilt, right/left	20°
Gallbladder support (option)	130 mm
Longitudinal movement (option)	230 mm

jörg&sohn

GmbH

Lindenstraße 24

DE-96482 Ahorn/Schorkendorf

Germany

Tel.: +49 (0) 95 61 / 82 99 0

Fax: +49 (0) 95 61 / 82 99 37

E-mail: info@joergundsohn.de

Internet: www.joergundsohn.de